



# PAT ACR Phantom Analysis

## Understanding the Report

### Report Header


The header of the report is at the top of the first page. Example below (Figure 2)

Facility: ATLANTA VAMC	Phantom: ACR	Contrast: 2.4:1
Scanner Model: SIEMENS Biograph64	Scan: 12/15/2011	Time Per Bed: 3.5min.
Reconstruction: PSF+TOF 2i21s Gauss3.00		

Figure 2. The Report Header

This Section reads the facility name, scanner make and model, reconstruction, scan date, and time per bed position from the DICOM tags. It also reports the actual contrast in the phantom based upon the reported activity injected into the phantom, spheres, and the phantom volume.

### Quantitative Assessment



At left is a PAT generated ACR analysis sheet (Figure 3). It is formatted identically to the standard ACR data sheet that is supposed to be filled in by the user.

Even though 20 quantitative measurements are required by ACR to be performed by the user, reviewers only look at three pieces of data.

Contrast				
	Hot Vial 8mm	Hot Vial 12mm	Hot Vial 16mm	Hot Vial 25mm
<b>max. SUV</b>	1.34	2.56	2.83	3.11 ✖

Scatter/Attenuation				
	Background	Bone	Air	Water
<b>mean SUV</b>	1.05 ✔	0.15	0.22	0.18
<b>min. SUV</b>	n/a	0.15	0.08	0.09

Ratio Calculations: max. vial SUV to mean background SUV				
	8mm/bkdg	12mm/bkdg	16mm/bkdg	25mm/bkdg
<b>Contrast</b>	1.27	2.43	2.68	2.95

Ratio Calculations: max. vial SUV to max. 25 mm vial			
	8mm/25mm	12mm/25mm	16mm/25mm
<b>Contrast</b>	0.43	0.82	0.91 ✔

Ratio Calculations: min. air or water to min. bone		
	air/bone	water/bone
<b>ratio</b>	1.98	2.04

Figure 3. ACR PET /CT phantom analysis

First, the “Hot Vial 25 cm”  $SUV_{max}$  must be between 1.8 and 2.8. If the value falls within this range, a green check mark is automatically generated. A red X is automatically placed in that box in Figure 3 because it is outside of ACR’s acceptance bracket. Interestingly, this data was taken on a well calibrated new SiPM array PET/CT scanner using the manufacturer’s most advanced clinical reconstruction. As this reconstruction appears to overcorrect for partial volume (as many of these advanced reconstructions do), it results in an out of range value.

Second, the  $SUV_{mean}$  of the background (equal to 1.0 if the fill was performed correctly) must be between 0.85 and 1.15. This is an overly generous acceptance for calibration, as every other accrediting body only allows a  $\pm 10\%$  calibration error, and most PET scanners can easily maintain a  $\pm 5\%$  accuracy. With this in mind, the PAT software generates a green check mark if the

# PAT ACR Phantom Analysis

calibration is within the  $\pm 10\%$  window, but a yellow warning icon is displayed if the  $SUV_{mean}$  is between 10-15% discrepant from the expected calibration. If the calibration is off by more than 10%, the scanner should be recalibrated.

The third quantitative measure is the ratio of the 16 mm cylinder  $SUV_{max}$  to the  $SUV_{max}$  of the 25 mm cylinder. This ratio must be greater than 0.7. If it meets this criterium, a green check mark is generated.

None of the other 17 quantitative measurements have been officially reviewed by ACR reviewers for more than a decade. As PAT stores all of this data electronically, it affords an opportunity to collect this data in an organized fashion, and to perhaps determine more empirically derived criteria for acceptable performance in the future.

### Image Uniformity/Calibration

A crude assessment of image uniformity is performed using the small uniform region in the center of the phantom. This is of limited utility given that only a handful of slices with uniform background are available for this assessment. This region is also used for the calibration assessment. An example of PAT output of Calibration and Uniformity is illustrated in Figure 4.




Figure 4. Calibration and Uniformity output.

## PAT ACR Phantom Analysis

### System Resolution/High contrast Resolution

Image resolution is to be determined from the smallest resolution pattern resolved in a 1 cm thick transaxial slice. A representative 1 cm thick slice is provided.




Figure 5. A representative 1 cm thick slice from which the smallest resolved rod pattern can be identified.

### Image Quality/Lesion Detectability




Figure 6. 1 cm thick slice used for low-contrast lesion detectability.

The low contrast lesion detectability test requires The 12 mm diameter cylinder (ROI 6 in Figure 6) must be visually discernable by the reviewer to meet ACR acceptance criteria.

# PAT ACR Phantom Analysis

## DICOM and Fill Information

Relevant DICOM header and fill information is displayed in fourth box (Figure 7). This is included to provide a simple means to check the fill and reconstruction information.

Name	Value
Institution	ATLANTA VAMC
Phantom	ACR
Series Description	PET WB 4min TrueX+TOF
Scan Date	12/15/2011
Scan Time	15:48:21
Assay Time	14:47:43
Patient Dose	10
Background Volume	5660.0g
Background Activity	0.83
Uptake Time	84.35
Cylinder Solution Volume	1000.0g
Cylinder Activity	0.35
Actual Contrast	2.40:1
Minutes per Bed	3.50
Voxel Dimensions	4.07x4.07x4.00mm
Matrix Dimensions	200x200x90
Scanner Make and Model	SIEMENS Biograph64
Reconstruction Method	PSF+TOF 2i21s
Reconstruction Parameters	
Reconstruction Filter	XYZ Gauss3.00

Figure 7. DICOM header and ACR Phantom fill information